
TMK 1 MODULAR REVISION STEM C E M E N T E D

TMK 1 MODULAR REVISION STEM

The key in

orthopaedic

sciences

**SURGICAL TECHNIQUE
CERRAHI TEKNİK**
IMPLANTS & INSTRUMENTS
İMLANTLAR VE EL ALETLERİ

CONTENTS İÇİNDEKİLER

Introduction
Giriş

2

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Introduction

Giriş

EN Used for hip arthroplasty in cases of hip fractures when collapse of the neck structure has occurred and revision arthroplasty where there is significant loss of bone stock in the region of the calcar. The calcar series offers 38-48 and 58 mm resection levels to compensate for bone deficiencies in the proximal/medial portion of the femur. In addition to the three resection levels, the stems are offered in three different diameters 8-10-12 mm. Diameter 8 mm component has standard 140 mm stem length, other stem diameters 10-12 mm offered with 150 mm stem length. These components are specifically designed to conform and match with the patient's anatomy. The keel, located below the medial collar provides the additional rotational stability that has been reduced due to medial bone loss. Three dimensional Interlock grit blast all over on the stem surface provides increased fixation at the prosthesis/cement interface compared with that of traditional smooth surface implants. The angle of the neck is 135° for all size.

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Material

Malzeme

EN It is manufactured from casting CoCr alloy according to ASTM F 75 - ISO 5832-4. This alloy is characterized by demonstrating excellent biocompatibility together with a high resistance of corrosion and superiormechnical strenght.

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Please Note : This document is intended as a guide for the surgeon only. There are multiple techniques for the insertion of TMK-1 Revision Stem(Cemented) and, as with any surgical procedure, a surgeon should be thoroughly trained and beware that this procedure is appropriate for the patient before proceeding.

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
EN The proximal femoral cuts can be fashioned by utilizing the trial prosthesis as a guide. The TMK-1 prosthesis comes in 23, 33, 43 and 53mm calcar resection levels. In general, the 23 and 33 mm calcar is utilized in situations where there has been little, if any, bone destruction in the proximal femur and a platform can be maintained at approximately 1cm above the lesser trochanter. The 43 mm design is utilized for bony defects down to the lower level of the lesser trochanter and the 53 mm prosthesis for defects below the lesser trochanter.

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
EN The TMK-1 prosthesis comes in 140mm and 150mm lengths. For 140mm and 150mm stems, tapered reamers are used to prepare the diaphyseal region. Medullary Awl reaming starts with using Ø8 Awl reamer (Ref: 10610321299) and if necessary, continues with Ø10 and Ø12 Awl reamers(Ref: 10610321301-1321). Reaming process finishes at least 2.0cm below the final implant level.

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EN When dealing with the compromised femur, it is recommended not to broaching or taking extra care during broaching. If the proximal femur is not entirely compromised, a rasp can be used to further shape the proximal envelope. Broaching is starts with the smallest rasp and if necessary, it is carried out with bigger sizes until the corresponding size of the last reamer is reached.

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EN The trial prosthesis should be very stable when it is driven in place. At this stage if necessary to provide desired anteversion Forked Femoral Inserter (Ref:10607020101) may be used. Once the trial is in position, mark where the keel comes into contact with the medial calcar. After marking the location for the keel, remove the trial stem and prepare a groove or recess in this same location on the medial bone using a Canal Chisel (Ref: 10605230011). After the groove has been fashioned, reinsert the trial stem, making sure it seats firmly on the horizontal femoral platform.

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EN With the trial calcar prosthesis in place, a trial reduction can be performed to evaluate proper leg length and lateral offset.

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**Attention/Dikkat**

EN A trial reduction cannot be performed with the broach. This will not provide accurate leg length and offset measurements.

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EN A distal cement restrictor is placed in the canal to allow a 2.0cm cement column below the tip of the stem. Cement is injected into the canal in a retrograde fashion and pressurized.

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! Attention/Dikkat

EN When using the TMK-1 implants, one must downsize the implant to allow for an adequate cement mantle. Example: Ream and broach to 12mm. Implant an 10mm TMK-1 stem.

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


EN The stem can be inserted with the Impactor Pim (Ref:10607031311). Place medial keel into medial calcar slot. Extraneous cement is then removed. When cement hardening has been achieved, the appropriate modular head determined during trial reduction can be impacted onto the stem (Ref: 10607020003).

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Forked Femoral Inserter
Çatallı Protez Çakıcı

10607020003

Femoral Head Impactor
Protez Çakıcı

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
Notes/Notlar

The products being manufactured by TIPSAN A.S. has been certified from SZUTEST (Strojirensk•
Zkušební Ústav) for EC certificate **(CE 1015)** according to Annex II. 3 and II.4 of Council Directive
93/42/EEC concerning medical devices

Kemalpasa Cd. 7404/1 Sk. No: 3 Pinarbasi - IZMIR - TURKEY
Tel: +90 232 479 56 54 Fax: +90 232 479 58 27
marketing@tipsan.com.tr www.tipsan.com.tr